REF CDI510H

CDI™ Shunt Sensor with heparin treatment
Третиран с хепарин паралелен сензор CDI™
Senzor za šant CDI™ obrađen heparinom
Paralelní senzor CDI™ s heparinovou úpravou
CDI™ Shuntsensor med heparinbehandling
CDI™ shuntsensor met heparinebehandeling
CDI™ hepariinraviga šundiandur
CDI™ -sunttianturi hepariinikäsittelyllä
Capteur pour shunt CDI™ traité à l'héparine
CDI™-Shunt-Sensor mit Heparinbehandlung
Αισθητήρας επικοινωνίας CDI™ με επεξεργασία ηπαρίνης
CDI™ heparinnal bevont shuntérzékeltő
CDI™ Shunt Sensor with heparin treatment
Sensore di derivazione CDI™ con trattamento eparina
CDI™ гепаринмен өңдеу мүмкіндігі бар шунттау датчигі
CDI™ šunta sensors ar heparīna pārklājumu
Heparinu apdorotas šunto jutiklis CDI™
Sensur tax-Shunt tas-CDI™ bi trattament bi-heparina
CDI™ shuntsensor med heparinbehandling
Czujnik bocznikowy CDI™ po obróbce heparyną
Sensor em derivação CDI™ com tratamento de heparina
Senzor de şunt CDI™ cu tratament cu heparină
Датчик для шунтирования CDI™ с обработкой гепарином
Senzor za šant CDI™ s heparinskim tretmanom
Paralelný senzor CDI™ impregnovaný heparínom
Šantni senzor CDI™ s heparinsko obdelavo
Sensor de derivación CDI™ con tratamiento de heparina
CDI™ shuntsensor med heparinbehandling
Heparin tedavisi ile CDI™ Şant Sensörü
Датчик-шунт CDI™ з обробкою гепарином

Table of Contents

ENGLISH.....	3
БЪЛГАРСКИ.....	11
HRVATSKI	21
ČESKY	29
DANSK	37
NEDERLANDS	45
EESTI KEEL	54
SUOMI	62
FRANÇAIS	70
DEUTSCH	79
ΕΛΛΗΝΙΚΑ	88
MAGYAR.....	98
IRISH.....	107
ITALIANO	115
ҚАЗАҚ ТІЛІ.....	124
LATVISKI	133
LIETUVIŲ KALBA.....	141
MALTI.....	149
NORSK	157
POLSKI	165
PORTUGUÊS.....	174
ROMÂNĂ.....	183
РУССКИЙ	192
SRPSKI	202
SLOVENČINA	210
SLOVENŠČINA.....	218
ESPAÑOL	226
SVENSKA	235
TÜRKÇE	244
УКРАЇНСЬКОЮ МОВОЮ	252

Symbols Glossary

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Symbol	Title	Description	Source
	Consult instructions for use or consult electronic instructions for use	Indicates the need for the user to consult the instructions for use	ISO 15223-1-5.4.3
	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified	ISO 15223-1-5.1.5
	Prescription only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician or other licensed practitioner	21 CFR 801.109
	Use – by date	Indicates the date after which the medical device is not to be used	ISO 15223-1-5.1.4
	Date of manufacture	Indicates the date when the medical device was manufactured	ISO 15223-1-5.1.3
	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified	ISO 15223-1-5.1.6
	Sterilized using irradiation	Indicates a medical device that has been sterilized using irradiation	ISO 15223-1-5.2.4
	Non-pyrogenic	Indicates a medical device that is non-pyrogenic	ISO 15223-1-5.6.3
	Fluid path	Indicates the presence of a fluid path	ISO 15223-1-5.6.2
	Manufacturer	Indicates the medical device manufacturer, as defined in EU Directives 90/385/EEC, 93/42/EEC, and 98/79/EC	ISO 15223-1-5.1.1

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PRODUCT DESCRIPTION

The CDI™ Shunt Sensor contains the K⁺, PO₂, PCO₂, and pH fluorescent microsensors, as well as the thermistor contact site for temperature measurement.

The CDI Shunt Sensor contains porcine derived heparin. The heparin treatment provides a non-leaching heparinized surface. It is applied by a process which covalently binds heparin to the blood contacting surfaces of this device.

The heparin treatment improves the thromboresistance of the non-biological materials. Surfaces with heparin treatment have demonstrated stability in the presence of blood.

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Note: The CDI Shunt Sensor can be placed in the shunt/purge line at any time during priming or bypass, as long as there is fluid in the circuit. You must stop the fluid flow upstream of the sensor before insertion to avoid loss of fluid.

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Warning: Do not remove the sensor filter/sparger assembly from the CDI Shunt Sensor until you are ready to connect it to the circuit. The microsenors on the CDI Shunt Sensor must be kept moist. Exposure to room air for more than a few minutes can damage the sensor.

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Warning:

- Prime solutions containing acetate ions such as Isolyte-S, Normosol-R, or Plasmalyte-A can cause damage to the PCO₂ sensor. If the pH channel reads less than 7.00 after the sensor is placed in the circuit, you should either recirculate the prime solution using a CO₂-free sweep gas or add sufficient buffer to raise the pH of the prime above 7.00. Exposure to acetate-containing prime solutions below pH 7.00 for longer than a few minutes can cause significant PCO₂ inaccuracy.
- Make sure all luer lock connections are securely tightened before priming the shunt/purge line. Connections that are not secure can result in a leak.
- The presence of air bubbles in CDI Shunt Sensors can affect the results. Intermittent bubbles, once removed from the CDI Shunt Sensor, will not affect its long-term accuracy.

4. The CDI Blood Parameter Monitoring System is now ready for use. Upon selecting the OPERATE mode, the monitor will begin reporting blood parameter values.

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Symbols Glossary

The following symbols may appear in the labeling, marking, or display of the Terumo Cardiovascular Systems (TCVS) CDI Shunt Sensor. These symbols are in accordance with the internationally harmonized standards.

Symbol	Title	Description	Source
	Consult instructions for use or consult electronic instructions for use	Indicates the need for the user to consult the instructions for use	ISO 15223-1-5.4.3
	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified	ISO 15223-1-5.1.5
Rx Only	Prescription only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician or other licensed practitioner	21 CFR 801.109
	Use – by date	Indicates the date after which the medical device is not to be used	ISO 15223-1-5.1.4
	Date of manufacture	Indicates the date when the medical device was manufactured	ISO 15223-1-5.1.3
	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified	ISO 15223-1-5.1.6
	Sterilized using irradiation	Indicates a medical device that has been sterilized using irradiation	ISO 15223-1-5.2.4
	Non-pyrogenic	Indicates a medical device that is non-pyrogenic	ISO 15223-1-5.6.3
	Fluid path	Indicates the presence of a fluid path	ISO 15223-1-5.6.2
	Manufacturer	Indicates the medical device manufacturer, as defined in EU Directives 90/385/EEC, 93/42/EEC, and 98/79/EC	ISO 15223-1-5.1.1

Symbol	Title	Description	Source
	Authorized representative in the European Community	Indicates the authorized representative in the European Community	ISO 15223-1-5.1.2
	Do not re-sterilize	Indicates a medical device that is not to be re-sterilized	ISO 15223-1-5.2.6
	Do not use if package is damaged and consult instructions for use	Indicates that a medical device that should not be used if the package has been damaged or opened and that the user should consult the instructions for use for additional information	ISO 15223-1-5.2.8
	Fragile, handle with care	Indicates a medical device that can be broken or damaged if not handled carefully	ISO 15223-1-5.3.1
	Keep dry	Indicates a medical device that needs to be protected from moisture	ISO 15223-1-5.3.4
	Do not re-use	Indicates a medical device that is intended for one single use only	ISO 15223-1-5.4.2
	Single sterile barrier system	Indicates a single sterile barrier system	ISO 11607-1-Annex E ISO 15223-1-5.2.11
	Medical device	Indicates that the device is a medical device	ISO 15223-1-5.7.7
	Unique device identifier	Indicates a carrier that contains unique device identifier information	ISO 15223-1-5.7.10
	Contents	Indicates the number of devices within the packaging.	NA
	Importer	Indicates the entity importing the medical device into the locale	ISO 15223-1-5.1.8
	Temperature limit	Indicates the temperature limits to which the medical device can be safely exposed.	ISO 15223-1-5.3.7

Symbol	Title	Description	Source
	Contains biological material of animal origin	Indicates a medical device that contains biological tissue, cells, or their derivatives, of animal origin	ISO 15223-1-5.4.8
	Contains a medical substance	Indicates a medical device that contains or incorporates a medical substance	ISO 15223-1-5.4.7

PRODUCT DESCRIPTION

The CDI™ Shunt Sensor contains the K⁺, PO₂, PCO₂, and pH fluorescent microsenors, as well as the thermistor contact site for temperature measurement.

The CDI Shunt Sensor contains porcine derived heparin. The heparin treatment provides a non-leaching heparinized surface. It is applied by a process which covalently binds heparin to the blood contacting surfaces of this device. The heparin treatment improves the thromboresistance of the non-biological materials. Surfaces with heparin treatment have demonstrated stability in the presence of blood.


Figure 1: The CDI Shunt Sensor with heparin treatment, Model CDI510H

INTENDED USE / PURPOSE

The disposable CDI Shunt Sensor with heparin treatment is intended for use with the CDI Monitoring Systems during cardiopulmonary bypass procedures when continuous blood gas, pH, and potassium monitoring is desired for up to 6 hours.

INTENDED USER

The CDI Shunt Sensor is intended for use by Healthcare Professionals trained to perform or assist in cardiopulmonary bypass procedures (e.g. Certified Clinical Perfusionist).

TARGETED PATIENT GROUPS

Patients undergoing cardiopulmonary bypass procedures that meet the minimum recommended device flow requirement of 35ml/min.

CLINICAL BENEFITS

Continuous in-line monitoring of blood gas parameters enables continuous patient monitoring to support perfusionist decisions during the CPB procedure and may indirectly improve patient outcomes.

INDICATIONS FOR USE

The CDI Shunt Sensor is indicated for use in cardiopulmonary bypass procedures to provide continuous blood gas, pH and potassium monitoring.

CONTRAINDICATIONS

The CDI510H is not intended for use in situations in which there is no flow of blood through the extracorporeal circuit. A minimum blood flow of 35 ml/min is recommended for optimal measurement performance of the shunt sensor. Restoration of minimum blood flow through the shunt line after an interruption will restore optimal performance of the system.

TRAINING REQUIREMENTS

This device is intended to be used by a medical professional working independently under the direction of a physician. The CDI Blood Parameter Monitoring System Operator's Manual contains complete instructions for assembly and use of the Shunt Sensors with the CDI Blood Parameter Monitoring System. It must be read in its entirety before attempting to assemble and use the system. To arrange additional training, please contact your local Terumo Cardiovascular Systems representative or call 1-800-521-2818 and ask about CDI Blood Monitoring Training.

WARNINGS

- **The CDI Shunt Sensors are sterile, heparin-coated, non-toxic, non-pyrogenic, single use devices, and for use in cardiopulmonary bypass procedures for up to 6 hours. Use beyond 6 hours may result in inaccuracies of the blood parameters and/or affect functional integrity of the device.**
- **The CDI Shunt Sensor is heparin treated and should not be used with heparin sensitive patients. Devices with heparin treated surfaces may cause an adverse reaction.**
- **Store CDI Shunt Sensors between 0°C (32°F) and 35°C (94°F). Freezing of the CDI Shunt Sensor, or storage at temperatures outside the stated range, can result in inaccurate performance. Check the freeze indicator on the box/carton before using the enclosed Shunt Sensors. If the clear indicator bulb has turned to violet, exposure to freezing temperatures has occurred and Shunt Sensors should not be used.**

- This device is sterilized using gamma irradiation and is intended for single use only. Do not reuse. Do not resterilize. Do not reprocess. Reprocessing may compromise the sterility, biocompatibility and functional integrity of the device.
- When administering intravascular dyes and novel pharmacological agents, or when dyshemoglobins or elevated bilirubin levels are present, independent external blood gas and blood chemistry analysis are required for accurate determination of all measured parameters needed to guide therapeutic decisions, as they can potentially cause inaccuracies in displayed value. Refer to the Operator's Manual for the CDI Blood Parameter Monitoring System for further specification of agents known to interfere with the system, the values affected and additional instructions.
- Maintain adequate levels of anticoagulation during extracorporeal circulation by monitoring activated clotting time (ACT) or other appropriate measurement. Use of a heparin treated device does not substitute for adequate anticoagulation levels.
- Verify the accuracy of displayed values with another source (i.e. laboratory or point of care blood gas analyzer) before initiating treatment.
- The CDI Shunt Sensor requires a minimum of 35 ml/min. Flows below the minimum may result in a slower time response. To maintain the minimum blood flow through the sensor, keep total blood flow in the shunt bypass line above 1.5 L/min for 1/2 inch tubing, 0.6 L/min for 3/8 inch tubing, and 0.2 L/min for 1/4 inch tubing. Restoration of blood flow above the minimum through the CDI Shunt Sensor will restore performance of the system.
- Use aseptic technique when inserting the shunt sensor (and shunt bypass line if used) into the extracorporeal circuit to ensure that the blood contact surfaces remain sterile.
- For all shunt bypass line applications: Use an arterial filter distal to the shunt bypass line when utilizing the shunt bypass line on the arterial side of the circuit. This protects against the introduction of air into the blood circulation.
- Exposure of the shunt sensor to prime solutions and/or blood with pH less than 7.0 or greater than 7.8 pH units or sodium measurement less than 120 or greater than 160 mEq/L can interfere in the accurate measurement of potassium.
- Failure to perform a 2 point tonometered gas calibration of the shunt sensor and potassium sensor may inhibit the system from achieving accuracy limits.

PRECAUTIONS

- Federal (USA) law restricts this device to sale by or on the order of a physician or other licensed practitioner.
- Do not use the CDI Shunt Sensor after the date printed on the package label. Use beyond the date may result in inaccurate performance.
- Do not use a CDI Shunt Sensor if the foil pouch it is packaged in has been damaged. A damaged foil pouch can result in inaccurate performance.

- Do not open the pouch until the Shunt Sensor is to be used. Exposure of the Shunt Sensor to room air for periods of greater than 24 hours may result in inaccurate calibration.
- Failure to follow the instructions for use can cause the monitoring system to display inaccurate values. The accuracy of the results is dependent upon the following:
 - Reading and understand the instructions for use
 - Proper setup and calibration of the system
 - Use of all available system features
 - Periodic comparison to a laboratory reference sample
 - Maintaining adequate blood flow

When a displayed value is significantly different from expectations based on the clinical situation, verify its accuracy by independent means before initiating treatment.
- The CDI Shunt Sensor contains Germall II in the calibration fluid. A potential byproduct of Germall II may be Formaldehyde. Exposure may cause adverse reactions in patients with Formaldehyde sensitivity.
- Do not connect a shunt sensor to an unprimed circuit. Prolonged “dry” exposure can damage the shunt sensors.

METHOD OF OPERATION

Shunt Sensor Placement

You need the following items: the Shunt Sensor(s) you are installing and access to the bypass circuit.

Note: CDI Shunt Sensor can be used for either arterial or venous application.

Warning: Once calibrated, do not remove and replace the shunt sensor from the cable-head prior to use. Removal and replacement could affect the measurement accuracy of the system.

Caution:

- Care should be taken to prevent BPM cable-heads from dropping onto any hard surface or otherwise receiving severe shock. If the cable-head had a shunt sensor installed and calibrated, you should replace the shunt sensor and repeat the calibration.
- The CDI Shunt Sensor should be placed distal to a purge line one-way valve to avoid possible back flow of air.
- Terumo recommends that you not connect the CDI Shunt Sensor directly onto another rigid plastic piece, such as a manifold, without extra support for the shunt sensor. Connecting the unsupported CDI Shunt Sensor directly to another rigid plastic piece may make the unsupported piece or the CDI Shunt Sensor susceptible to breakage. Terumo recommends having a length of pliable tubing between a supported CDI Shunt Sensor and any other rigid plastic part.

Note: The CDI Shunt Sensor can be placed in the shunt/purge line at any time during priming or bypass, as long as there is fluid in the circuit.

You must stop the fluid flow upstream of the sensor before insertion to avoid loss of fluid.

Note: The sensor is bidirectional. Blood can flow through it in either direction.

Note: If you are placing the CDI Shunt Sensor in a sampling line, place it on the inlet side of the sample port to avoid intermittent interruption of blood parameter data during medication administration.

Note: The end of the CDI Shunt Sensor that has the large blue luer cap is a male luer connector. The end of the CDI Shunt Sensor attached to the filter/sparger assembly is a female luer connector.

Follow these steps to install the CDI Shunt Sensor into a shunt/purge line:

1. Using sterile technique, remove the top luer cap (white) from the CDI Shunt Sensor and attach one end of the “shunt/purge” line to the top of the CDI Shunt Sensor. Make sure the large blue top venting luer is completely tightened.
2. Remove the sensor filter/sparger assembly from the bottom of the sensor. Attach the other end of the shunt/purge tubing circuit to the sensor.

Warning: Do not remove the sensor filter/sparger assembly from the CDI Shunt Sensor until you are ready to connect it to the circuit. The microsensors on the CDI Shunt Sensor must be kept moist. Exposure to room air for more than a few minutes can damage the sensor.


Figure 2: Installation into the circuit

3. Prime and debubble the shunt/purge line, inspecting the line and CDI Shunt Sensor for bubbles. Air bubbles may clear more easily from the CDI Shunt Sensor if it is in an upright position. Intermittent bubbles, once removed from the CDI Shunt Sensor, will not affect its long-term accuracy. Inspect the primed line for leaks; do not use the sensor if any leaks are detected.

Warning:

- Prime solutions containing acetate ions such as Isolyte-S, Normosol-R, or Plasmalyte-A can cause damage to the PCO₂ sensor. If the pH channel reads less than 7.00 after the sensor is placed in the circuit, you should either recirculate the prime solution using a CO₂-free sweep gas or add sufficient buffer to raise the pH of the prime above 7.00. Exposure to acetate-containing prime solutions below pH 7.00 for longer than a few minutes can cause significant PCO₂ inaccuracy.
- Make sure all luer lock connections are securely tightened before priming the shunt/purge line. Connections that are not secure can result in a leak.

- **The presence of air bubbles in CDI Shunt Sensors can affect the results. Intermittent bubbles, once removed from the CDI Shunt Sensor, will not affect its long-term accuracy.**
4. The CDI Blood Parameter Monitoring System is now ready for use. Upon selecting the OPERATE mode, the monitor will begin reporting blood parameter values.

Shunt Sensor Replacement

Shunt Sensor Replacement During Cardiopulmonary Bypass: Terumo does not recommend replacing the shunt sensor during cardiopulmonary bypass. In the event that the CDI Shunt Sensor(s) or the CDI monitor fails, the case should be completed using routine laboratory blood gas analysis.

DISPOSAL

After use, dispose of product and packaging in accordance with hospital, administrative and/or local, state, federal and international laws and regulations.

STORAGE

Temperature: 0°C to 35°C (32°F to 95°F)

SPECIFICATIONS

Priming volume: 1.2 mL

Product Complaints

Any serious incident that has occurred in relation to the device should be reported to Terumo and the competent authority of the Member State in which the user and/or patient is established.

Potential Undesirable Side Effects

Delay to procedure, adverse bodily reaction, infection, hemolysis, embolism, blood loss, unable to meet performance claims.

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**TERUMO** and CDI are trademarks of Terumo and may be registered.
© Copyright Terumo Cardiovascular Systems Corporation 2005-2008
AU SPONSOR TERUMO AUSTRALIA PTY LTD.